

Anti-MonoMethyl-Histone H3 (Lys9) (10F5) Mouse antibody



Catalog Number: 168154

Product name

WB

Anti-MonoMethyl-Histone H3 (Lys9) (10F5) Mouse antibody

Dilutions

WB: 1/1000-3000

Specificity

Human, Mouse, Rat

Antibody description

MonoMethyl-Histone H3 (Lys9) (10F5) Mouse monoclonal antibody

Preparation

Antigen: Synthetic Peptide of Histone H3 (Mono Methyl Lys9)

Formulation

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality


Monoclonal

Ig Type


Mouse IgG1

Applications

Validations



Western blot analysis of 1) Rat Testis Tissue, 2) Raw264.7, diluted at 1:1000. cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).



Western blot detection of Mono-Methyl-Histone H3 in MCF7, C2C12, C6, 3T3 cell lysates using Mono-Methyl-Histone H3 (Lys9) (10F5) Mouse mAb (1:1000 diluted). Predicted band size: 15KDa. Observed band size: 15KDa.