

Anti-Flotillin 1 (6H9) Mouse antibody



Catalog Number: 168168

Product name

Anti-Flotillin 1 (6H9) Mouse antibody

Specificity

Mouse, Rat

Antibody description

Flotillin 1 (6H9) Mouse monoclonal antibody

Preparation

Antigen: Synthetic Peptide of Flotillin-1

Formulation

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB


Dilutions

WB: 1/1000-2000

Validations



Western blot analysis of 1) Mouse Brain Tissue, 2) Rat Brain tissue, diluted at 1:2000.



Western blot detection of Flotillin-1 in Rat Brain, Mouse Brain cell lysates using Flotillin-1 (6H9) Mouse mAb(1:1000 diluted). Predicted band size: 49KDa. Observed band size: 49KDa.