

# Anti-Nucleophosmin (2D9) Mouse antibody



Catalog Number: 168179

## Product name

Anti-Nucleophosmin (2D9) Mouse antibody

## Ig Type

Mouse IgG1

## Specificity

Human, Mouse

## Applications

WB, ChIP

## Antibody description

Nucleophosmin (2D9) Mouse monoclonal antibody

## Dilutions

WB: 1/1000

## Preparation

Antigen: Purified recombinant full length of human NPM1 protein expressed in E.coli.

;ICC: 1/200

## Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Validations



Western blot detection of NPM1 in HeLa, Jurkat and 3T3 cell lysates using NPM1 mouse mAb (1:1000 diluted). Predicted band size: 33KDa. Observed band size: 38KDa.