

Anti-TriMethyl-Histone H3 (Lys79) (9G4) Mouse antibody



Catalog Number: 168187

Product name

WB

Anti-TriMethyl-Histone H3 (Lys79) (9G4) Mouse antibody

Dilutions

WB: 1/500-2000

Specificity

Human, Mouse, Rat

Antibody description

TriMethyl-Histone H3 (Lys79) (9G4) Mouse monoclonal antibody

Preparation

Antigen: Synthetic Peptide of Histone H3 (Tri Methyl Lys79)

Formulation

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality


Monoclonal

Ig Type


Mouse IgG1

Applications

Validations



Western blot detection of Tri-Methyl-Histone H3 in MCF7,C6,Hela cell lysates using Tri-Methyl-Histone H3 (Lys79) (9G4) Mouse mAb(1:1000 diluted).Predicted band size:15KDa.Observed



Western blot analysis of HeLa, diluted at 1:500. cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003,Inventbiotech,MN,USA).