

# Anti-Histone H3 (Zebrafish Specific) (2C4) Mouse antibody



Catalog Number: 168188

## Product name

WB

Anti-Histone H3 (Zebrafish Specific) (2C4) Mouse antibody

## Dilutions

WB: 1/1000-3000

## Specificity

Zebrafish

## Antibody description

Histone H3 (Zebrafish Specific) (2C4) Mouse monoclonal antibody

## Preparation

Antigen: Synthetic Peptide of Histone H3

## Formulation

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality


Monoclonal

## Ig Type


Mouse IgG1

## Applications

## Validations



Western blot detection of Histone H3 in C6,3T3,HeLa cell lysates using Histone H3 (Zebrafish Specific) (2C4) Mouse mAb(1:10000 diluted).Predicted band size:15KDa.Observed



band size:15KDa.

Western blot analysis of Zebrafish skeletal muscle, diluted at 1:2000. cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003,Inventbiotech,MN,USA).