

# Anti-Peroxiredoxin 1/PAG (6A6) Mouse antibody



Catalog Number: 168206

## Product name

Anti-Peroxiredoxin 1/PAG (6A6) Mouse antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Peroxiredoxin 1/PAG (6A6) Mouse monoclonal antibody

## Preparation

Antigen: Recombinant Protein of Peroxiredoxin-1

## Formulation

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

WB, IHC-P, ICC/IF


## Dilutions

WB: 1/1000-3000


IF: 1/100-200

IHC: 1/50-300


## Validations



Western blot detection of Prdx 1 in T47D,3T3,COS7,C6 and HeLa cell lysates using Prdx 1 mouse mAb (1:3000 diluted). Predicted band size:21kDa. Observed band size:21kDa.



Western blot detection of Peroxiredoxin 1/Prdx1 in Jurkat,CHO-K1,C6,HeLa cell lysates using Peroxiredoxin 1/Prdx1 (6A6) Mouse mAb(1:1000 diluted). Predicted band size:21KDa. Observed band size:21KDa.




IF analysis of HeLa with antibody (Left) and DAPI (Right) diluted at 1:100.

# Anti-Peroxiredoxin 1/PAG (6A6) Mouse antibody



Catalog Number: 168206



Western blot analysis of 1) MCF7, 2) Rat Kidney Tissue, 3) Mouse Brain Tissue, diluted at 1:2000.