

# Anti-Placental alkaline phosphatase (PLAP) Mouse antibody



Catalog Number: 168207

## Product name

Anti-Placental alkaline phosphatase (PLAP) Mouse antibody

## Specificity

Human

## Antibody description

Placental alkaline phosphatase (PLAP) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant Placental alkaline phosphatase (PLAP) protein fragments expressed in E.coli.

## Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

WB

## Dilutions

WB: 1/1000

## Validations

