

Anti-HDAC3 (3G3) Mouse antibody



Catalog Number: 168208

Product name

Anti-HDAC3 (3G3) Mouse antibody

Specificity

Human, Mouse, Rat

Antibody description

HDAC3 (3G3) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human HDAC3 protein fragments expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB


Dilutions

WB: 1/1000


Validations



Western blot detection of HDAC3 in THP-1, LNCAP, MCF7, PC12, 293, Jurkat, 3T3 and HeLa cell lysates using HDAC3 mouse mAb (1:1000 diluted). Predicted band size: 49 kDa. Observed band size: 49 kDa.



Western blot analysis of extracts from HDAC3 knock down 293T, wild-type (WT), K562, Jurkat and HeLa cell lysates using HDAC3 Mouse mAb (upper, 1:500 diluted) or Beta Actin (200068-8F10) Mouse mAb (lower, 1:10000 diluted). Predicted band size: 49 kDa. Observed band size: 49 kDa.



Western blot detection of HDAC3 in Jurkat and HeLa cell lysates using HDAC3 mouse mAb

Anti-HDAC3 (3G3) Mouse antibody



Catalog Number: 168208

(1:1000 diluted).Predicted band

size:49KDa.Observed band size:49KDa.