

# Anti-BRG1 (4E11) Mouse antibody



Catalog Number: 168211

## Product name

Anti-BRG1 (4E11) Mouse antibody

## Specificity

Human

## Antibody description

BRG1 (4E11) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human BRG1 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

WB, IP


## Dilutions

WB: 1/1000

## Validations



Western blot detection of SBRG1 in HeLa NE and K562 cell lysates using BRG1 mouse mAb (1:1000 diluted). Predicted band size: 220KDa. Observed



band size: 220KDa.

Immunoprecipitation analysis of HeLa cells lysates. IP: BRG1 mouse mAb(168211); WB: BRG1 mouse mAb(201025).