

Anti-HAUSP/USP7 (3E4) Mouse antibody



Catalog Number: 168213

Product name

Anti-HAUSP/USP7 (3E4) Mouse antibody

Specificity

Human

Antibody description

HAUSP/USP7 (3E4) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human HAUSP / USP7 protein fragments expressed in E.coli.

Formulation

Ascites containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB

Dilutions

WB: 1/200

Validations



Western blot detection of HAUSP / USP7 in HeLa and K562 cell lysates using HAUSP / USP7 mouse mAb (1:200 diluted), with super ECL. Predicted band size: 135KDa. Observed band size: 135KDa.