

Anti-S100A4 (3D10) Mouse antibody



Catalog Number: 168218

Product name

WB, IP

Anti-S100A4 (3D10) Mouse antibody

Dilutions

Specificity

WB: 1/500

Human

Antibody description

S100A4 (3D10) Mouse monoclonal antibody

Preparation

Antigen: Recombinant human S100A4 protein.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality


Monoclonal

Ig Type


Mouse IgG1

Applications

Validations



Immunoprecipitation analysis of HeLa cell lysates



using S100A4 mouse mAb.

Western blot detection of S100A4 in HeLa and A549 cell lysates using S100A4 mouse mAb(dilution 1:500).Predicted band size:12kDa.Observed band size:12kDa.