

Anti-Liver FABP (8G6) Mouse antibody



Catalog Number: 168220

Product name

ELISA

Anti-Liver FABP (8G6) Mouse antibody

Validations

Specificity

Human

Antibody description

Liver FABP (8G6) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human L-FABP protein fragments expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300.

Storage

Store at 4°C short term.Avoid freeze / thaw cycle.


Clonality

Monoclonal


Ig Type

Mouse IgG1

Applications



Observed L-FABP levels in several Human urine samples and Human serum samples at different dilution.



Standard Curve for L-FABP: Capture Antibody Mouse mAb (8G6-D5) to L-FABP and Detector Antibody Mouse mAb(2D6-F6) to L-FABP.