

# Anti-TNF alpha (5F3) Mouse antibody



Catalog Number: 168223

## Product name

Anti-TNF alpha (5F3) Mouse antibody

## Specificity

Human

## Antibody description

TNF alpha (5F3) Mouse monoclonal antibody

## Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications

ELISA

## Validations



Standard Curve for TNF-a: Capture Antibody Mouse mAb 168223 to TNF-a at 4ug/ml and Anti-hTNF Human antibody for detecting.