

Anti-TNF alpha (5F3) Mouse antibody



Catalog Number: 168223

Product name

Anti-TNF alpha (5F3) Mouse antibody

Ig Type

IgG1

Specificity

Human

Applications

ELISA

Antibody description

TNF alpha (5F3) Mouse monoclonal antibody

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Standard Curve for TNF-a: Capture Antibody Mouse mAb 168223 to TNF-a at 4ug/ml and Anti-hTNF Human antibody for detecting.