

# Anti-SMYD4 (7C8) Mouse antibody



Catalog Number: 168226

## Product name

Anti-SMYD4 (7C8) Mouse antibody

## Specificity

Human

## Antibody description

SMYD4 (7C8) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human SMYD4 protein fragments expressed in E.coli.

## Formulation

Ascites containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

WB


## Dilutions

WB: 1/200;Other 1/100

## Validations



Western blot detection of SMYD4 in 40ug HeLa, 40ug MCF7 and 40ug MDA-MB-468 cell lysates using SMYD4 mouse mAb (1:100 diluted). Predicted band size: 89KDa. Observed



band size: 89KDa.

Western blot detection of SMYD4 in HeLa, MCF7 and A431 cell lysates using SMYD4 mouse mAb (1:200 diluted). Predicted band size: 89KDa. Observed band size: 89KDa.