

# Anti-TRIM21 (5B9) Mouse antibody



Catalog Number: 168227

## Product name

Anti-TRIM21 (5B9) Mouse antibody

## Specificity

Human

## Antibody description

TRIM21 (5B9) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human TRIM21 protein fragments expressed in E.coli.

## Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of TRIM21 in Embryonic, MCF7, 293T and HeLa cell lysates using TRIM21 mouse mAb (dilution 1:1000). Predicted band size: 50 kDa. Observed band size: 50 kDa.