

# Anti-CARM1 (2B9) Mouse antibody



Catalog Number: 168237

## Product name

WB, IP

Anti-CARM1 (2B9) Mouse antibody

## Dilutions

## Specificity

WB: 1/200-1/500

Human, Mouse

## Antibody description

CARM1 (2B9) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human PRMT4/CARM1 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine(pH 7.4 150 mM NaCl)with 0.02% sodium azide 0.1mg/mlBSA and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality


Monoclonal

## Ig Type


IgG1

## Applications

## Validations



Immunoprecipitation analysis of HeLa cell lysates using PRMT4/CARM1 mouse mAb.



Western blot detection of PRMT4/CARM1 in HeLa,A431 and K562 cell lysates using PRMT4/CARM1 mouse mAb (1:200-1:500 diluted).Predicted band size:63KDa.Observed band size:63KDa.