

Anti-XRCC1 (3F6) Mouse antibody



Catalog Number: 168242

Product name

Anti-XRCC1 (3F6) Mouse antibody

Ig Type

IgG1

Specificity

Human

Applications

WB

Antibody description

XRCC1 (3F6) Mouse monoclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: Purified recombinant human XRCC1 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Western blot detection of XRCC1 in Raji, K562, HeLa, Jurkat and THP-1 cell lysates using XRCC1 mouse mAb (1:1000 diluted). Predicted band size: 82KDa. Observed band size: 82KDa.