

Anti-Hexokinase II (3B9) Mouse antibody



Catalog Number: 168247

Product name

Anti-Hexokinase II (3B9) Mouse antibody

Ig Type

Mouse IgG1

Specificity

Human

Applications

WB

Antibody description

Hexokinase II (3B9) Mouse monoclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: Purified recombinant human Hexokinase II protein fragments expressed in E.coli.

Validations

Formulation


PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal



Western blot detection of Hexokinase II in COS7, K562 and HeLa cell lysates using Hexokinase II mouse mAb (1:1000 diluted). Predicted band size: 102 kDa. Observed band size: 102 kDa.