

Anti-Wnt3a (6H11) Mouse antibody



Catalog Number: 168248

Product name

Anti-Wnt3a (6H11) Mouse antibody

Ig Type

Mouse IgG1

Specificity

Human

Applications

WB

Antibody description

Wnt3a (6H11) Mouse monoclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: Purified recombinant Wnt3a protein fragments expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Western blot detection of Wnt3a in MCF-7,ZR75-1,MDA-MB-468 and T47D cell lysates using Wnt3a mouse mAb (1:1000 diluted).Predicted band size: 42KDa.Observed band size:42KDa.