

# Anti-SIRT6 (6C9) Mouse antibody



Catalog Number: 168257

## Product name

Anti-SIRT6 (6C9) Mouse antibody

## Specificity

Human, Mouse, Rat, Monkey

## Antibody description

SIRT6 (6C9) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human SIRT6 protein expressed in E.coli.

## Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


WB, ICC/IF, IP

## Dilutions


WB: 1/500

;ICC: 1/100


## Validations



Immunofluorescent analysis of HeLa cells fixed by 4% paraformaldehyde and using SIRT6 mouse mAb (dilution 1:100).



Immunoprecipitation analysis of HeLa cell lysates using SIRT6 mouse mAb.



Western blot analysis of extracts from Jurkat, MCF7, HCT116, K562, COS7, C2C12, 3T3 and HeLa cell lysates using SIRT6 mouse mAb (1:500 diluted). Predicted band size: 42,36kDa. Observed band size: 42,36kDa.