

Anti-Fatty Acid Synthase (6B9) Mouse antibody



Catalog Number: 168259

Product name

Anti-Fatty Acid Synthase (6B9) Mouse antibody

Specificity

Human, Mouse, Rat, Monkey

Antibody description

Fatty Acid Synthase (6B9) Mouse monoclonal antibody

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


WB, IP, ICC/IF

Dilutions


WB: 1/1000

;ICC: 1/200


Validations



Immunoprecipitation analysis of CHO-K1 cell lysates using Fatty Acid Synthase mouse mAb.



Western blot detection of Fatty Acid Synthase in HeLa, C6, 3T3, CHO-K1 and MCF7 cell lysates using Fatty Acid Synthase mouse mAb (dilution 1:500). Predicted band size: 273 kDa. Observed band



size: 273 kDa.

Immunocytochemistry staining of HeLa cells fixed with 4% Paraformaldehyde and using anti-Fatty Acid Synthase mouse mAb (dilution 1:200).