

Anti-ALDH2 (8D4) Mouse antibody



Catalog Number: 168262

Product name

Anti-ALDH2 (8D4) Mouse antibody

Specificity

Human, Mouse, Rat, Monkey, Hamster

Antibody description

ALDH2 (8D4) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human ALDH2 protein fragments expressed in E.coli

Formulation

Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of ALDH2 in HeLa, 3T3, C6, COS7, CHO-K1 and A549 cell lysates using ALDH2 mouse mAb (1:1000 diluted). Predicted band size: 56 kDa. Observed band size: 56 kDa.