

# Anti-Lipocalin-2/NGAL (3C6) Mouse antibody



Catalog Number: 168265

## Product name

Monoclonal

Anti-Lipocalin-2/NGAL (3C6) Mouse antibody

## Ig Type

## Specificity

IgG2a

Human

## Applications

## Antibody description


ELISA

Lipocalin-2/NGAL (3C6) Mouse monoclonal antibody

## Validations

## Preparation

Antigen: Purified recombinant human NGAL protein expressed in E.coli.



## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300.

## Storage

Store at 4°C short term.Avoid freeze / thaw cycle.

Anti-FITC mouse mAb was used to capture FITC-labeled NGAL Mouse mAb (168265), HRP-labeled NGAL Mouse mAb (201179-2F6) was used for detecting.

## Clonality