

Anti-Anthrax Protective gen (5C1) Mouse antibody



Catalog Number: 168268

Product name

Anti-Anthrax Protective gen (5C1) Mouse antibody

Ig Type

Mouse IgG2b

Specificity

Bacillus anthracis

Applications

ELISA

Antibody description

Anthrax Protective gen (5C1) Mouse monoclonal antibody

Validations

Formulation


PBS(pH 7.4) containing with 0.03% Proclin300 and 25% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal



Standard Curve for Anthrax Protective Antigen: Capture Antibody Mouse mAb 168268 to Anthrax Protective Antigen at 4u03bcg/ml and other antibody for detecting.