

Anti-Visfatin (2E11) Mouse antibody



Catalog Number: 168269

Product name

Anti-Visfatin (2E11) Mouse antibody

Specificity

Human

Antibody description

Visfatin (2E11) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of PBEF1 (aa338-479) expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


WB

Dilutions

WB: 1/500 - 1/2000

ELISA: 1/10000

Validations



Western blot detection of PBEF/NAMPT in Jurkat, A549, MCF7, C6 and HeLa cell lysates using PBEF/NAMPT mouse mAb (1:5000 diluted). Predicted band size: 56 kDa. Observed band size: 52 kDa.