

Anti-IkB alpha (2E10) Mouse antibody



Catalog Number: 168272

Product name

WB

Anti-IkB alpha (2E10) Mouse antibody

Dilutions

Specificity

WB: 1/1000;Other 1/500

Human, Mouse

Antibody description

IkB alpha (2E10) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human IkB-alpha(N-terminus) protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality


Monoclonal

Ig Type


Mouse IgG

Applications

Validations



Western blot detection of IkB-alpha(N-terminus) in 3T3 cell lysate using IkB-alpha(N-terminus) mouse mAb (1:500 diluted). Predicted band size: 36KDa. Observed band size: 36KDa.



Western blot detection of IkB-alpha(N-terminus) in Jurkat and HeLa cell lysates using IkB-alpha(N-terminus) mouse mAb (1:1000 diluted). Predicted band size: 36KDa. Observed band size: 36KDa.