

# Anti-JNK1 (1A4) Mouse antibody



Catalog Number: 168281

## Product name

Anti-JNK1 (1A4) Mouse antibody

## Specificity

Human, Mouse, Rat

## Antibody description

JNK1 (1A4) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human JNK1 protein fragments expressed in E.coli.

## Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications


WB, ICC/IF

## Dilutions


WB: 1/1000

;ICC: 1/100


## Validations



Immunofluorescent analysis of 3T3 cells fixed by anhydrous methanol for 2 h at -20°C and using anti-JNK1 mouse mAb (dilution 1:100).



Western blot detection of JNK1 in C6 and 3T3 cell lysates using JNK1 mouse mAb (1:1000 diluted). Predicted band size: 46,54KDa. Observed



band size: 46KDa.

Western blot detection of JNK1 in CHO-K1 cell lysate (B) and CHO-K1 transfected by JNK1-fragment fusion protein (A) cell lysate using JNK1 mouse mAb (1:2000 diluted). Predicted band size: 46,54KDa. Observed band size: 46,54KDa.