

Anti-TBL1 (4H2) Mouse antibody



Catalog Number: 168288

Product name

Anti-TBL1 (4H2) Mouse antibody

Dilutions

WB: 1/1000

Specificity

Human

;ICC: 1/100


Antibody description

TBL1 (4H2) Mouse monoclonal antibody

Validations

Preparation

Antigen: Purified recombinant human TBL1 protein fragments expressed in E.coli.



Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Western blot detection of TBL1X in SW480 cell lysates using TBL1X mouse mAb (1:1000 diluted). Predicted band size:58KDa. Observed band size:58KDa.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality


Monoclonal

Ig Type

Mouse IgG2b

Applications

WB, ICC/IF



Immunocytochemistry staining of HeLa cells fixed with 4% Paraformaldehyde and using anti-TBL1X mouse mAb (dilution 1:100).