

Anti-RBBP5 (1C9) Mouse antibody



Catalog Number: 168290

Product name

Anti-RBBP5 (1C9) Mouse antibody

Specificity

Human, Mouse, Rat

Antibody description

RBBP5 (1C9) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human RBBP5 protein fragments expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


WB, ICC/IF

Dilutions


WB: 1/1000

;ICC: 1/100

Validations



Western blot detection of RBBP5 in C2C12, THP-1, PC12 and 3T3 cell lysates using RBBP5 mouse mAb (1:1000 diluted). Predicted band size: 70KDa. Observed band size: 70KDa.



Immunofluorescent analysis of HeLa cells fixed with 4% Paraformaldehyde and using anti-RBBP5 mouse mAb (dilution 1:100). DAPI was used to stain nucleus (blue).