

Anti-Vimentin (4F2E8) Mouse antibody



Catalog Number: 168297

Product name

Mouse IgG1

Anti-Vimentin (4F2E8) Mouse antibody

Applications

Specificity

WB, IHC-P

Human, Monkey

Dilutions

Antibody description

WB: 1/500 - 1/2000

Vimentin (4F2E8) Mouse monoclonal antibody

IHC: 1/200 - 1/1000

Preparation

ELISA: 1/10000

Antigen: Purified recombinant fragment of Vimentin (aa2-466) expressed in E. Coli.

Validations

Formulation

Ascitic fluid containing 0.03% sodium azide.


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of Vimentin in U2OS,293,COS7 and HeLa cell lysates using Vimentin mouse mAb (dilution 1:500).Predicted band size:54KDa.Observed band size:54KDa.