

# Anti-SMN/Gemin 1 (2F1) Mouse antibody



Catalog Number: 168301

## General Information

### Product name:

Anti-SMN/Gemin 1 (2F1) Mouse antibody

### Immunogen:

### Applications:

WB, IHC-P, FC, IP

### Specificity:

Human, Monkey

### Antibody description:

SMN/Gemin 1 (2F1) Mouse monoclonal antibody

## Properties

### Preparation:

Antigen: Purified recombinant fragment of human SMN1 expressed in E. Coli.

### Formulation:

Ascitic fluid containing 0.03% sodium azide.

### Storage:

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

### Clonality:

Monoclonal

**Ig Type:** Mouse IgG1

## Applications


WB: 1/500 - 1/2000

IHC: 1/200 - 1/1000


FC: 1/200 - 1/400

ELISA: 1/10000

## Validations



Immunoprecipitation analysis of HeLa cell lysates using Gemin 1 mouse mAb.



Western blot detection of Gemin 1 in COS7, Jurkat and HeLa cell lysates using Gemin 1 mouse mAb (dilution 1:1000). Predicted band size: 35 Kda. Observed band size: 39KDa.