

# Anti-C17orf53 (7E1) Mouse antibody



Catalog Number: 168326

## Product name

Mouse IgG2a

Anti-C17orf53 (7E1) Mouse antibody

## Applications

## Specificity

WB

Human

## Dilutions

## Antibody description

WB: 1/500 - 1/2000

C17orf53 (7E1) Mouse monoclonal antibody

ELISA: 1/10000


## Preparation

Antigen: Purified recombinant fragment of human C17ORF53 (AA: 282-527) expressed in E. Coli.

## Validations

## Formulation

Purified antibody in PBS with 0.05% sodium azide



## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of C17orf53 in Jurkat, A549, MCF7 and HeLa cell lysates using C17orf53 mouse mAb (1:1000 diluted). Predicted band size: 70KDa. Observed band size: 70KDa.

## Clonality

Monoclonal

## Ig Type