

Anti-NRF1 (9E6) Mouse antibody



Catalog Number: 168328

General Information

Product name:

Anti-NRF1 (9E6) Mouse antibody

Immunogen:

Applications:

WB

Specificity:

Human

Antibody description:

NRF1 (9E6) Mouse monoclonal antibody

Properties

Preparation:

Antigen: Purified recombinant human NRF1 protein fragments expressed in E.coli

Formulation:

Purified mouse monoclonal in PBS(pH 7.4)containing with 0.02% sodium azide 50% glycerol 0.1%BSA.

Storage:

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality:


Monoclonal

Ig Type: IgG2a

Applications

WB: 1/1000

Validations



Western blot detection of NRF1 in HeLa, K562, CEM and U2OS cell lysates using TCF11/NRF1 mouse mAb (1:1000 diluted). Predicted band size: 67 kDa. Observed band size: 72 kDa.