

Anti-ALDH10/ALDH3A2 (3F7) Mouse antibody



Catalog Number: 168329

Product name

Monoclonal

Anti-ALDH10/ALDH3A2 (3F7) Mouse antibody

Ig Type

Specificity

Mouse IgG2a

Human

Applications

Antibody description

WB

ALDH10/ALDH3A2 (3F7) Mouse monoclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: Purified recombinant human Aldehyde dehydrogenase 10 protein fragments expressed in E.coli.

Formulation


Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Validations



Western blot detection of Aldehyde dehydrogenase 10 in A549, U251 and HepG2 cell lysates using Aldehyde dehydrogenase 10 mouse mAb (1:1000 diluted). Predicted band size: 55KDa. Observed band size: 55KDa.