

Anti-S100 beta (5G3) Mouse antibody



Catalog Number: 168333

Product name

Anti-S100 beta (5G3) Mouse antibody

Specificity

Rat

Antibody description

S100 beta (5G3) Mouse monoclonal antibody

Preparation

Antigen: Human S100B synthetic peptide conjugated to KLH.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG2b


Applications

WB

Dilutions

WB: 1/500

Validations



Western blot analysis of extracts from Rat Brain cell lysates using S100B mouse mAb (1:500 diluted). Predicted band size: 10kDa. Observed band size: 10kDa.