

Anti-E Cadherin (3H11) Mouse antibody



Catalog Number: 168334

Product name

Mouse IgG1

Anti-E Cadherin (3H11) Mouse antibody

Applications

Specificity

WB

Human, Mouse, Monkey

Dilutions

Antibody description

WB: 1/500 - 1/2000

E Cadherin (3H11) Mouse monoclonal antibody

ELISA: 1/10000

Preparation

Antigen: Purified recombinant fragment of human CDH1 expressed in E. Coli.

Validations

Formulation

Ascitic fluid containing 0.03% sodium azide.


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of E-Cadherin in PC3, Caco2 and HCT116 cell lysates using E-Cadherin mouse mAb (dilution 1:500). Predicted band size: 135 kDa. Observed band size: 135 kDa.