

# Anti-LC3A/B (3E9) Mouse antibody



Catalog Number: 168336

## Product name

Anti-LC3A/B (3E9) Mouse antibody

## Ig Type

IgG2b

## Specificity

Human, Rat

## Applications

WB

## Antibody description

LC3A/B (3E9) Mouse monoclonal antibody

## Dilutions

WB: 1/1000

## Preparation

Antigen: Synthetic peptide corresponding to human LC3B protein

## Validations

## Formulation


PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal



Western blot detection of LC3A/B in Rat Cerebral Cortex, Rat Brain and Hela cell lysates using LC3A/B Mouse mAb (1:1000 diluted). Predicted band size: 15KDa. Observed band size: 14, 16KDa.