

# Anti-IKB beta (4B11) Mouse antibody



Catalog Number: 168338

## Product name

Mouse IgG1

Anti-IKB beta (4B11) Mouse antibody

## Applications

## Specificity

WB

Human

## Dilutions

## Antibody description

WB: 1/500 - 1/2000

IKB beta (4B11) Mouse monoclonal antibody

ELISA: 1/10000

## Preparation

Antigen: Purified recombinant fragment of human NFKBIB expressed in E. Coli.

## Validations

## Formulation

Ascitic fluid containing 0.03% sodium azide.


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot detection of IκBβ in Jurkat cell lysates using IκBβ mouse mAb (1:1000 diluted). Predicted band size: 48KDa. Observed band size: 48KDa.