

# Anti-Cyclin H (4E11) Mouse antibody



Catalog Number: 168339

## Product name

WB, IP

Anti-Cyclin H (4E11) Mouse antibody

## Dilutions

## Specificity

WB: 1/1000

Human

## Antibody description

Cyclin H (4E11) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human Cyclin H protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality


Monoclonal

## Ig Type


Mouse IgG2b

## Applications

## Validations



Immunoprecipitation analysis of HeLa cell lysates using Cyclin H mouse mAb.



Western blot detection of Cyclin H in Jurkat, K562 and A431 cell lysates using Cyclin H mouse mAb (1:1000 diluted). Predicted band size: 38KDa. Observed band size: 38KDa.