

# Anti-EPO (7E1) Mouse antibody



Catalog Number: 168341

## Product name

Anti-EPO (7E1) Mouse antibody

## Specificity

Human

## Antibody description

EPO (7E1) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human EPO expressed in E. Coli.

## Formulation

Ascitic fluid containing 0.03% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ICC/IF

## Dilutions

ICC: 1/200 - 1/1000

ELISA: 1/10000

## Validations



Indirect ELISA assay for anti-EPO mouse mAb. Antigen coating concentration: 8ug/ml.