

Anti-BMPR2 (1F12) Mouse antibody



Catalog Number: 168357

Product name

Mouse IgG1

Anti-BMPR2 (1F12) Mouse antibody

Applications

Specificity

WB, IHC-P

Human, Mouse, Rat, Monkey

Dilutions

Antibody description

WB: 1/500 - 1/2000

BMPR2 (1F12) Mouse monoclonal antibody

IHC: 1/200 - 1/1000

Preparation

ELISA: 1/10000

Antigen: Purified recombinant fragment of human BMPR2 expressed in E. Coli.

Validations

Formulation

Ascitic fluid containing 0.03% sodium azide.


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of BMPR2 in COS7, HeLa and 3T3 cell lysates using BMPR2 mouse mAb (dilution 1:500). Predicted band size: 115 Kda. Observed band size: 115KDa.