

Anti-PP-X (2F11) Mouse antibody



Catalog Number: 168359

Product name

Anti-PP-X (2F11) Mouse antibody

Specificity

Human

Antibody description

PP-X (2F11) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human Protein Phosphatase 4C protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


WB, IHC-P, ICC/IF

Dilutions


WB: 1/200

;ICC: 1/200


Validations



Immunocytochemistry of HeLa cells using anti-Protein Phosphatase 4C mouse mAb diluted 1:200.



IHC of paraffin-embedded human breast cancer using anti-Protein Phosphatase 4C mouse mAb



diluted 1/500-1/1000.

Western blot detection of Protein Phosphatase 4C in HeLa and Jurkat cell lysates using Protein Phosphatase 4C mouse mAb (1:200 diluted). Predicted band size: 34KDa. Observed band size: 34KDa.