

Anti-Dynein Light Chain Mouse antibody



Catalog Number: 168370

Product name

Anti-Dynein Light Chain Mouse antibody

Specificity

Human

Antibody description

Dynein Light Chain Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human DNAL4 (AA: 1-105) expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, FC


Dilutions

WB: 1/500 - 1/2000


FC: 1/200 - 1/400

ELISA: 1/10000


Validations



Western blot analysis using DNAL4 mAb against HEK293 (1) and DNAL4 (AA: 1-105)-hIgGfC transfected HEK293 (2) cell lysate.



Flow cytometric analysis of HEK293 cells using DNAL4 mouse mAb (green) and negative control



(purple).

Western blot analysis of extracts from C6 cell lysates using Dynein light chain mouse mAb (1:500 diluted). Predicted band size: 12KDa. Observed band size: 12KDa.