

Anti-BRAF (4E1) Mouse antibody



Catalog Number: 168376

Product name

Anti-BRAF (4E1) Mouse antibody

Ig Type

Mouse IgG2a

Specificity

Human, Mouse

Applications

WB

Antibody description

BRAF (4E1) Mouse monoclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: Purified recombinant human B Raf protein fragments expressed in E.coli.

Validations

Formulation


Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal



Western blot detection of B Raf in MDA-MB-468, T47D, MCF7, HCT116, C2C12, 3T3 and HeLa cell lysates using B Raf mouse mAb (1:1000 diluted). Predicted band size: 87 kDa. Observed band size: 87 kDa.