

Anti-STAT5a (5C10) Mouse antibody



Catalog Number: 168378

Product name

Anti-STAT5a (5C10) Mouse antibody

Specificity

Human, Mouse

Antibody description

STAT5a (5C10) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human Stat5a protein fragments expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2b


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of Stat5a in 3T3, A549, Molt-4, K562 and HeLa cell lysates using Stat5a Mouse mAb (1:1000 diluted). Predicted band size: 90KDa. Observed band size: 90KDa.