

# Anti-Bax (2B12) Mouse antibody



Catalog Number: 168380

## Product name

Anti-Bax (2B12) Mouse antibody

## Specificity

Human, Monkey

## Antibody description

Bax (2B12) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human Bax protein fragments expressed in E.coli

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

IgG2a


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of Bax in COS7,Raji,CEM and HeLa cell lysates using Bax mouse mAb (1:1000 diluted). Predicted band size:20KDa. Observed band size:20KDa.