

# Anti-SMYD1 (2E11) Mouse antibody



Catalog Number: 168385

## Product name

Anti-SMYD1 (2E11) Mouse antibody

## Specificity

Rat

## Antibody description

SMYD1 (2E11) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human SMYD1 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b


## Applications

WB

## Dilutions

WB: 1/100

## Validations



Western blot detection of SMYD1 in rat skeletal muscle lysates and using SMYD1 mouse mAb (1:100 diluted). Predicted band size: 56KDa. Observed band size: 56KDa.