

Anti-Human IgG (9G6) Mouse antibody



Catalog Number: 168389

Product name

Anti-Human IgG (9G6) Mouse antibody

Specificity

Human

Antibody description

Human IgG (9G6) Mouse monoclonal antibody

Preparation

Antigen: Human IgG was isolated from human sera and purified by chromatography.

Formulation

Purified antibody in PBS containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


WB

Dilutions

WB: 1/500 - 1/2000

ELISA: 1/10000

Validations



Western blot detection of human IgG using human IgG mouse mAb (1:1000 diluted). Predicted band size: 50KDa. Observed band size: 50KDa.