

# Anti-SWAP70 (2E9) Mouse antibody



Catalog Number: 168391

## Product name

Anti-SWAP70 (2E9) Mouse antibody

## Specificity

Human, Mouse, Rat, Monkey

## Antibody description

SWAP70 (2E9) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human SWAP70 protein fragments expressed in E.coli.

## Formulation

Ascites containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b


## Applications

WB, IP


## Dilutions

WB: 1/1000

## Validations



Immunoprecipitation analysis of HeLa cell lysate using SWAP70 mouse mAb.



Western blot detection of SWAP70 in COS7,293T,C6,3T3 and HeLa cell lysates and using SWAP70 mouse mAb (1:1000 diluted).Predicted band size: 70KDa.Observed band size: 70KDa.